










# Size and weight chart VINCENTevolution5



## VINCENTevolution5 SIZES


					
	XS	S	M	L	XL

WRIST OPTIONS WITH LAMINATION RING/PLATE						
<b>VINCENTwrist quicksnap</b> Lamination ring coaxial 6 	<b>Total length</b>	170 mm 6.69"	180 mm 7.08"	190 mm 7.48"	200 mm 7.87"	210 mm 8.26"
	<b>Total weight</b>	ca. 477 g	ca. 481 g	ca. 533 g	ca. 535 g	ca. 575 g
<b>VINCENTwrist quicksnap + flexion</b> Lamination ring coaxial 6 	<b>Total length</b>	189 mm 7.44"	199 mm 7.83"	209 mm 8.22"	219 mm 8.62"	229 mm 9.01"
	<b>Total weight</b>	ca. 536 g	ca. 540 g	ca. 592 g	ca. 594 g	ca. 634 g
<b>VINCENTwrist short</b> Transcarpal radial 4 	<b>Total length</b>	146 mm 5.74"	156 mm 6.14"	166 mm 6.53"	176 mm 6.92"	186 mm 7.32"
	<b>Total weight</b>	ca. 422 g	ca. 426 g	ca. 478 g	ca. 480 g	ca. 520 g
<b>VINCENTwrist short + flexion</b> Transcarpal radial 4 	<b>Total length</b>	165 mm 6.49"	175 mm 6.88"	185 mm 7.28"	195 mm 7.67"	205 mm 8.07"
	<b>Total weight</b>	ca. 481 g	ca. 485 g	ca. 537 g	ca. 539 g	ca. 579 g

# Size and weight chart VINCENTyoung3+



## VINCENTyoung3+



### WRIST OPTIONS WITH LAMINATION RING/PLATE

<p><b>VINCENTwrist quicksnap</b> <i>Lamination ring koaxial 6</i></p>	<b>Total length</b>	166 mm 6.53 "
	<b>Total weight</b>	approx. 391 g
<p><b>VINCENTwrist quicksnap + flexion</b> <i>Lamination ring koaxial 6</i></p>	<b>Total length</b>	185 mm 7.28 "
	<b>Total weight</b>	approx. 436 g
<p><b>VINCENTwrist short</b> <i>Transcarpal radial 4</i></p>	<b>Total length</b>	142 mm 5.59 "
	<b>Total weight</b>	approx. 347 g
<p><b>VINCENTwrist short + flexion</b> <i>Transcarpal radial 4</i></p>	<b>Total length</b>	161 mm 6.33 "
	<b>Total weight</b>	approx. 387 g