










Size and weight chart VINCENTevolution5



VINCENTevolution5 SIZES

					
	XS	S	M	L	XL

WRIST OPTIONS WITH LAMINATION RING/PLATE						
VINCENTwrist quicksnap Lamination ring coaxial 6 	Total length	170 mm 6.69"	180 mm 7.08"	190 mm 7.48"	200 mm 7.87"	210 mm 8.26"
	Total weight	ca. 477 g	ca. 481 g	ca. 533 g	ca. 535 g	ca. 575 g
VINCENTwrist quicksnap + flexion Lamination ring coaxial 6 	Total length	189 mm 7.44"	199 mm 7.83"	209 mm 8.22"	219 mm 8.62"	229 mm 9.01"
	Total weight	ca. 536 g	ca. 540 g	ca. 592 g	ca. 594 g	ca. 634 g
VINCENTwrist short Transcarpal radial 4 	Total length	146 mm 5.74"	156 mm 6.14"	166 mm 6.53"	176 mm 6.92"	186 mm 7.32"
	Total weight	ca. 422 g	ca. 426 g	ca. 478 g	ca. 480 g	ca. 520 g
VINCENTwrist short + flexion Transcarpal radial 4 	Total length	165 mm 6.49"	175 mm 6.88"	185 mm 7.28"	195 mm 7.67"	205 mm 8.07"
	Total weight	ca. 481 g	ca. 485 g	ca. 537 g	ca. 539 g	ca. 579 g